

Product datasheet for **MP207687**

Lpcat1 Mouse qPCR Primer Pair (NM_145376)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	210992
Forward Sequence:	CCTGTATGCCAGCAATGTGAGG
Reverse Sequence:	TCAGCAGGCAAGCGAAGTTGTC
ACCN:	<u>BC066809</u> , <u>NM_145376</u> , <u>NM_145376.1</u> , <u>NM_145376.2</u> , <u>NM_145376.3</u> , <u>NM_145376.4</u> , <u>NM_145376.5</u> , <u>BC005662</u> , <u>NM_145376.6</u>
UniProt ID:	<u>Q3TFD2</u>
Synonyms:	2900035H07Rik; Ayt12; BB137372; BC005662; C87117; LPCAT; rd11
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)