

## Product datasheet for **MP207698**

### Lamp3 Mouse qPCR Primer Pair (NM\_177356)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	239739
Forward Sequence:	GTCCTAAATGGAAGCAGGCTCTG
Reverse Sequence:	GGAGTCACATTTCCAGATGCAT
ACCN:	<u>NM_177356</u> , <u>NM_177356.1</u> , <u>NM_177356.2</u> , <u>NM_177356.3</u> , <u>BC160231</u> , <u>BY067461</u> , <u>BY388271</u>
UniProt ID:	<u>Q7TST5</u>
Synonyms:	1200002D17Rik; Cd208; DC-LAMP; DCLAMP; LAMP; TSC403
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)