

Product datasheet for **MP207721**

Gbp6 Mouse qPCR Primer Pair (NM_194336)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	100702
Forward Sequence:	TCGAGAGTTCCATCTTGCAGGC
Reverse Sequence:	AGCAGCTCTTGCTCCTTCTCTG
ACCN:	<u>BC007143</u> , <u>BC057969</u> , <u>BC115768</u> , <u>NM_194336</u> , <u>NM_194336.2</u> , <u>BC057969.1</u> , <u>BB657433</u> , <u>BC052867</u> , <u>BK005759</u> , <u>BQ562803</u> , <u>BX525828</u> , <u>BY159737</u>
Synonyms:	AI595338; Mpa2l
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)