

Product datasheet for **MP207747**

Magt1 Mouse qPCR Primer Pair (NM_025952)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67075
Forward Sequence:	CGGTGGATTGCAGACAGAACTG
Reverse Sequence:	CTTCGCAGATACACAAGTCCACC
ACCN:	<u>BC003881, NM_025952, NM_025952.1, NM_025952.2, NM_025952.3, NM_025952.4, BC023242</u>
UniProt ID:	<u>Q9CQY5</u>
Synonyms:	2410001C15Rik; 2610529C04Rik; 2810482I07Rik; IAG2; IAP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)