

Product datasheet for **MP207753**

Mfsd11 Mouse qPCR Primer Pair (NM_178620)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	69900
Forward Sequence:	GTCGCCATTGTTGGACCTCAGA
Reverse Sequence:	GACAGATGCTGTGTAGAAGGACC
ACCN:	BC088988 , NM_178620 , NM_178620.1 , NM_178620.2 , NM_178620.3 , BC016209 , BC079583 , BC100547 , BC137970 , BC137971 , BC145068 , NM_178620.4
UniProt ID:	Q8BJ51
Synonyms:	2600014M03Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)