

## Product datasheet for **MP207768**

### **Naa35 Mouse qPCR Primer Pair (NM\_030153)**

#### **Product data:**

|                   |                                                                                                                                                                                                                                                                                                                                               |
|-------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:     | qPCR Primer Pairs                                                                                                                                                                                                                                                                                                                             |
| Gene ID:          | 78689                                                                                                                                                                                                                                                                                                                                         |
| Forward Sequence: | CTGCTCACAGTGCTCATAGCCT                                                                                                                                                                                                                                                                                                                        |
| Reverse Sequence: | TGGTAGTGCCATTCTGAGCCTG                                                                                                                                                                                                                                                                                                                        |
| ACCN:             | <a href="#">BC056435</a> , <a href="#">NM_030153</a> , <a href="#">NM_030153.1</a> , <a href="#">NM_030153.2</a> , <a href="#">BC010475</a> , <a href="#">BC027201</a> , <a href="#">BC043701</a> , <a href="#">BY100294</a>                                                                                                                  |
| UniProt ID:       | <a href="#">Q6PHQ8</a>                                                                                                                                                                                                                                                                                                                        |
| Synonyms:         | A330021G12Rik; A330027C19Rik; AI158944; C030004C14Rik; Mak10                                                                                                                                                                                                                                                                                  |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.                                                                                                                                   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.                                                                                                                                                                                                                                                                                                                               |
| Stability:        | The primer mix is stable for one year from date of shipping.                                                                                                                                                                                                                                                                                  |



[View online »](#)