

## Product datasheet for **MP207769**

### Mkrn1 Mouse qPCR Primer Pair (NM\_018810)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54484
Forward Sequence:	TCAGAGGACTGGGTGAATGCCA
Reverse Sequence:	TTGTCTCCACCGTAGTTGGCTC
ACCN:	<u><a href="#">BC003329</a></u> , <u><a href="#">NM_018810</a></u> , <u><a href="#">NM_018810.1</a></u> , <u><a href="#">NM_018810.2</a></u> , <u><a href="#">BM933910</a></u> , <u><a href="#">BP761527</a></u>
Synonyms:	RFP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)