

## Product datasheet for **MP207823**

### Mgat5 Mouse qPCR Primer Pair (NM\_145128)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	107895
Forward Sequence:	CGCGAATGGATGTGCTTTCCTG
Reverse Sequence:	GACTTCTGCGTAAGGATGCTGG
ACCN:	<a href="#">BC125533</a> , <a href="#">BC125535</a> , <a href="#">NM_145128</a> , <a href="#">NM_145128.2</a> , <a href="#">NM_145128.3</a> , <a href="#">BC030320</a> , <a href="#">BC053924</a> , <a href="#">BC094225</a> , <a href="#">BY752687</a>
UniProt ID:	<a href="#">Q8R4G6</a>
Synonyms:	4930471A21Rik; 5330407H02Rik; AI480971; GlcNAc-TV
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)