

Product datasheet for **MP207837**

Map7d2 Mouse qPCR Primer Pair (NM_001081124)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	78283
Forward Sequence:	TGGAGGAAAGGCTGAGCACAGT
Reverse Sequence:	CCTCTTGTTCTTCTGTAGCCG
ACCN:	NM_001081124 , NM_001081124.1 , NM_001081124.2 , BE952779
UniProt ID:	A2AG50
Synonyms:	1600028E09Rik; 2900002G04Rik; 5330432J06Rik; C030036K04Rik; Mtap7d2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

