

Product datasheet for **MP207886**

Mmp14 Mouse qPCR Primer Pair (NM_008608)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	17387
Forward Sequence:	GGATGGACACAGAGAACTTCGTG
Reverse Sequence:	CGAGAGGTAGTTCTGGGTTGAG
ACCN:	<u>BC076638</u> , <u>NM_008608</u> , <u>NM_008608.1</u> , <u>NM_008608.2</u> , <u>NM_008608.3</u> , <u>NM_008608.4</u> , <u>BY700899</u>
UniProt ID:	<u>P53690</u>
Synonyms:	AI325305; MMP-X1; MT-MMP-1; MT1-MM; MT1-MMP; sa; sabe
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)