

Product datasheet for **MP207920**

Mkks Mouse qPCR Primer Pair (NM_021527)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	59030
Forward Sequence:	CACTGTCTGCACCTTGCTTCTC
Reverse Sequence:	CACCTCCAATAAAACCCACGG
ACCN:	BC117836 , NM_021527 , NM_021527.1 , NM_021527.2 , BC024359 , BC058174 , BC080765 , BC100336 , BC117835 , BY735964
UniProt ID:	Q9JL70
Synonyms:	1300013E18Rik; AI463362; AI957237; Bbs; Bbs6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)