

Product datasheet for **MP207925**

Med13l Mouse qPCR Primer Pair (NM_172424)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	76199
Forward Sequence:	CACGGAGTTTAGGATGGAAGTGG
Reverse Sequence:	AAGGCTGGAAGTCGGCACTTT
ACCN:	<u>NM_172424</u> , <u>NM_172424.1</u> , <u>NM_172424.2</u> , <u>NM_172424.3</u> , <u>NM_172424.4</u> , <u>BC028781</u> , <u>BC040283</u> , <u>BC049985</u> , <u>BC052320</u> , <u>BC072653</u> , <u>BC156485</u> , <u>BC172732</u> , <u>BQ748188</u>
UniProt ID:	<u>Q6JPI3</u>
Synonyms:	2210413I17Rik; 6330591G05Rik; 9030618F05Rik; mKIAA1025; Thrap2; Trap240L
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)