

Product datasheet for **MP207933**

Med23 Mouse qPCR Primer Pair (NM_027347)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	70208
Forward Sequence:	CACTACCGAGTTCAGCTCCTGA
Reverse Sequence:	CTGTGATGAGCCTCAGTGCAGT
ACCN:	<u>BC050916</u> , <u>NM_027347</u> , <u>NM_027347.1</u> , <u>NM_027347.2</u> , <u>NM_027347.3</u> , <u>BC005508</u> , <u>BC040398</u> , <u>BQ885176</u>
UniProt ID:	<u>Q80YQ2</u>
Synonyms:	130kDa; 3000002A17Rik; Crsp3; ESTM7; mKIAA1216; Sur2; X83317
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)