

Product datasheet for **MP207934**

Med24 Mouse qPCR Primer Pair (NM_011869)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23989
Forward Sequence:	CCTGTGGGACAAAGCTGACCAG
Reverse Sequence:	AGGCTGGCAAAGTTGACCTCAG
ACCN:	<u>BC005409</u> , <u>NM_011869</u> , <u>NM_011869.1</u> , <u>NM_011869.2</u> , <u>BC078459</u> , <u>BE950589</u> , <u>BU612260</u> , <u>BY034665</u> , <u>NM_011869.3</u>
UniProt ID:	<u>Q99K74</u>
Synonyms:	100kD; 911GSE; AU040102; AW547152; D11Ertd307; D11Ertd307e; DRIP; DRIP100; Gse; Gse2; Pp; Pparb2; Pparbp2; R7552; R75526; Thr; Thrap4; Trap; Trap100
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)