

Product datasheet for **MP207935**

Med26 Mouse qPCR Primer Pair (NM_027485)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	70625
Forward Sequence:	CCATCACCAAAGAGGCACTGGA
Reverse Sequence:	GTAGTGCCACTTCGTTCTGGTG
ACCN:	BC054737 , NM_027485 , NM_027485.1 , NM_027485.2 , NM_027485.3 , NM_027485.4 , BC024555 , BY059601
UniProt ID:	Q7TN02
Synonyms:	5730493L18Rik; AI414941; AW495270; Crsp7
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)