

Product datasheet for **MP207944**

Med12 Mouse qPCR Primer Pair (NM_021521)

Product data:

| | |
|-------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 59024 |
| Forward Sequence: | GCTCCTAGAAGACCTGATTCGC |
| Reverse Sequence: | AGGATTCCGGCAAGTAAGCCTGG |
| ACCN: | <u>NM_021521</u> , <u>NM_021521.1</u> , <u>NM_021521.2</u> , <u>BC057085</u> , <u>BC057119</u> , <u>BC090654</u> , <u>BC110696</u> , <u>BC145447</u> |
| UniProt ID: | <u>A2AGH6</u> |
| Synonyms: | 230kDa; Mopa; OPA-1; Tnrc11; Trap230 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)