

Product datasheet for **MP207962**

Meis2 Mouse qPCR Primer Pair (NM_010825)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	17536
Forward Sequence:	AGACCACGATGACGCAACCTCA
Reverse Sequence:	TGCCAGGTGAAGCTACGCTGT
ACCN:	BC017375 , NM_010825 , NM_010825.1 , NM_010825.2 , NM_010825.3 , BB612458 , BQ443762 , BQ572308 , BY012436 , BY290986
UniProt ID:	P97367
Synonyms:	A430109D20Rik; Mei; Mrg; Mrg1; Str; Stra10
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)