

Product datasheet for **MP207992**

Pwwp3b Mouse qPCR Primer Pair (NM_175541)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	245631
Forward Sequence:	GGACTCTGCTTCCATACTGCTG
Reverse Sequence:	CCTTTGTGTCCATCTCACTTGGC
ACCN:	BC120584 , NM_175541 , NM_175541.1 , NM_175541.2 , NM_175541.3 , BC053077 , BC063766 , BC068178 , BC096534 , BC137665 , BY124815 , BY713797 , BY729100
UniProt ID:	Q4VA55
Synonyms:	9430086K21Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)