

Product datasheet for **MP207997**

MIph Mouse qPCR Primer Pair (NM_053015)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	171531
Forward Sequence:	AGGATGCTCACGCTGAACTGGA
Reverse Sequence:	GAGGTGTCTACATCAGCCAGGT
ACCN:	<u>NM_053015</u> , <u>NM_053015.1</u> , <u>NM_053015.2</u> , <u>NM_053015.3</u> , <u>BC150804</u> , <u>BC006900</u>
UniProt ID:	<u>Q91V27</u>
Synonyms:	2210418F23Rik; 5031433I09Rik; AW228792; D1Wsu84e; I(1)-3Rk; I1Rk3; In; Slac-2a
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)