

Product datasheet for **MP208007**

Mboat7 Mouse qPCR Primer Pair (NM_029934)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	77582
Forward Sequence:	TGCTATGGCACAGACTTCTGCG
Reverse Sequence:	ACTCAACAGCATGGTCCAGGCA
ACCN:	<u>BC023417</u> , <u>BC118958</u> , <u>NM_029934</u> , <u>NM_029934.1</u> , <u>NM_029934.2</u> , <u>NM_029934.3</u>
UniProt ID:	<u>Q8CHK3</u>
Synonyms:	5730589L02Rik; BB1; Leng4; Lpiat; Lpiat1; mBB1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)