

Product datasheet for **MP208008**

Mme Mouse qPCR Primer Pair (NM_008604)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	17380
Forward Sequence:	CAGCCGAAACTACAAGGAGTCC
Reverse Sequence:	CATAAAGCCTCCCCACAGCATTCC
ACCN:	<u>BC034092</u> , <u>NM_008604</u> , <u>NM_008604.1</u> , <u>NM_008604.2</u> , <u>NM_008604.3</u> , <u>NM_008604.4</u> , <u>BC066840</u> , <u>BY145991</u>
UniProt ID:	<u>Q61391</u>
Synonyms:	6030454K05Rik; C85356; CALLA; CD10; NEP; SFE
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)