

Product datasheet for **MP208016**

Mpp7 Mouse qPCR Primer Pair (NM_001081287)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	75739
Forward Sequence:	AGAGTGAGCAGTATGACACGGC
Reverse Sequence:	TCGCTTCAGCTCATTTAGCCCG
ACCN:	<u>BC117550</u> , <u>BC118058</u> , <u>NM_001081287</u> , <u>NM_001081287.1</u> , <u>NM_001081287.2</u> , <u>BC049662</u>
UniProt ID:	<u>Q8BVD5</u>
Synonyms:	1110068J02Rik; 2810038M04Rik; 5430426E14Rik; AI415104; Gm955
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)