

Product datasheet for **MP208023**

Marchf5 Mouse qPCR Primer Pair (NM_027314)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	69104
Forward Sequence:	AGCCGTGACTTATGGAGCAGTG
Reverse Sequence:	TATCAGCATGACAGGAATAGTAGG
ACCN:	BC064752 , BC107215 , BC132455 , BC132457 , NM_027314 , NM_027314.1 , NM_027314.2 , NM_027314.3 , BC107214 , BY070559
UniProt ID:	Q3KNM2
Synonyms:	1810015H18Rik; 2310008I22Rik; 2700055A20Rik; 5730499H23Rik; E130202O05Rik; MARCH-V
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)