

Product datasheet for **MP208024**

Marchf6 Mouse qPCR Primer Pair (NM_172606)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	223455
Forward Sequence:	GTTCCCTCTCATTTGTGGCTGG
Reverse Sequence:	CCAGTGCAGAAACATCGTGGTC
ACCN:	BC059190 , NM_172606 , NM_172606.1 , NM_172606.2 , BC027068 , BC031721 , BC037454 , BC039776 , BC048816 , BU516606
UniProt ID:	Q6ZQ89
Synonyms:	3830408G03; F830029L24Rik; mKIAA0597
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)