

Product datasheet for **MP208035**

Ms4a2 Mouse qPCR Primer Pair (NM_013516)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	14126
Forward Sequence:	CCTCAATCTGACCAATAACTTCGC
Reverse Sequence:	GCCAGGATGGTGAGAAACAGCA
ACCN:	<u>NM_013516</u> , <u>NM_013516.1</u> , <u>NM_013516.2</u> , <u>BC139324</u> , <u>BB847820</u> , <u>BC139323</u> , <u>BY585879</u>
UniProt ID:	<u>P20490</u>
Synonyms:	Fce1b; Fcer1; Fcer1b; fcERI; FcR; FcRB; Fcrbe; Fcrlbeta; Ms4a1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)