

Product datasheet for **MP208056**

Mesd Mouse qPCR Primer Pair (NM_023403)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67943
Forward Sequence:	GCGATTACAACGATGCCGACATG
Reverse Sequence:	CGATAGGTGCTGAGGGTCTCTT
ACCN:	BC014742 , NM_023403 , NM_023403.1 , NM_023403.2 , NM_023403.3 , BC058612 , BU054434 , NM_023403.4
UniProt ID:	Q9ERE7
Synonyms:	2210015O11Rik; AW537813; mesd; mKIAA0081; msd
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)