

Product datasheet for **MP208065**

Met Mouse qPCR Primer Pair (NM_008591)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	17295
Forward Sequence:	GTTCTGCTTGGCAACGAGAGCT
Reverse Sequence:	GGAGAATGCACTGTATTGCGTCG
ACCN:	<u>NM_008591</u> , <u>NM_008591.1</u> , <u>NM_008591.2</u> , <u>BC117969</u> , <u>BC117970</u> , <u>BC128481</u> , <u>BC128482</u> , <u>BC156287</u> , <u>BU514228</u>
Synonyms:	AI838057; c-Met; HGF; HGFR; Par4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)