

Product datasheet for **MP208097**

Metap1 Mouse qPCR Primer Pair (NM_175224)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	75624
Forward Sequence:	GTGGATGAAGGAGCACGGAAC
Reverse Sequence:	CCATTTGCTTGGCGTGTTCTG
ACCN:	<u>BC096469, NM_175224, NM_175224.1, NM_175224.2, NM_175224.3, NM_175224.4, BC027110</u>
UniProt ID:	<u>Q8BP48</u>
Synonyms:	1700029C17Rik; AW047992; mKIAA0094
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)