

Product datasheet for **MP208107**

Mccc1 Mouse qPCR Primer Pair (NM_023644)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	72039
Forward Sequence:	AGAGCTGCGAAGGCTGTCAAGT
Reverse Sequence:	GGATGTTCCAATTGTAGCCTGG
ACCN:	<u>BC021382</u> , <u>NM_023644</u> , <u>NM_023644.1</u> , <u>NM_023644.2</u> , <u>NM_023644.3</u> , <u>NM_023644.4</u> , <u>BY139695</u>
UniProt ID:	<u>Q99MR8</u>
Synonyms:	1810045E08Rik; 2310058B18Rik; Mcca; MCCalpha; R75106
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)