

Product datasheet for **MP208110**

Mthfd1l Mouse qPCR Primer Pair (NM_172308)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	270685
Forward Sequence:	AGGAGGTGGATATGCTCAGGTC
Reverse Sequence:	GCCAGCAAGTTATTAGCGGCAG
ACCN:	BC030437 , NM_172308 , NM_172308.1 , NM_172308.2 , NM_172308.3 , NM_172308.4 , BC049936 , BC058265
UniProt ID:	Q3V3R1
Synonyms:	2410004L15Rik; AI647056; Fthfsdc1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)