

Product datasheet for **MP208143**

Mcph1 Mouse qPCR Primer Pair (NM_173189)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	244329
Forward Sequence:	CCACTGGATTTCTGAGGAGCCT
Reverse Sequence:	GAGGGTTCCTTGGTATTGCTGG
ACCN:	<u>BC056924</u> , <u>NM_173189</u> , <u>NM_173189.1</u> , <u>NM_173189.2</u> , <u>NM_173189.3</u>
UniProt ID:	<u>Q7TT79</u>
Synonyms:	5430437K10Rik; BRIT1; D030046N04Rik; MCT
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)