

Product datasheet for **MP208155**

Morc3 Mouse qPCR Primer Pair (NM_001045529)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	338467
Forward Sequence:	AATGGCAACGGGATGACAGCAG
Reverse Sequence:	CGCATAGAACCTGACTTGAAGCC
ACCN:	<u>BC098072</u> , <u>NM_001045529</u> , <u>NM_001045529.1</u> , <u>NM_001045529.2</u> , <u>NM_001045529.3</u> , <u>BC026506</u> , <u>BC029284</u> , <u>BC145705</u>
UniProt ID:	<u>F7BJB9</u>
Synonyms:	1110051N18Rik; AI452146; BF318192; D16Jhu32e; NXP2; Zcwcc3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)