

## Product datasheet for **MP208158**

### Mgst2 Mouse qPCR Primer Pair (NM\_174995)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	211666
Forward Sequence:	TCGTAATGCTGTGGATGGCTGG
Reverse Sequence:	TTTCTCAGCGGCTTCGGCATAG
ACCN:	<a href="#">BC132234</a> , <a href="#">BC132260</a> , <a href="#">NM_174995</a> , <a href="#">NM_174995.1</a> , <a href="#">NM_174995.2</a> , <a href="#">NM_174995.3</a> , <a href="#">BC145529</a> , <a href="#">BX639527</a>
UniProt ID:	<a href="#">A2RST1</a>
Synonyms:	GST2; MGST-II
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)