

Product datasheet for **MP208175**

Map6 Mouse qPCR Primer Pair (NM_010837)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	17760
Forward Sequence:	CTACCAAACCCGACGACAAGGA
Reverse Sequence:	GTGGCTCTTTGGTTACAGGACC
ACCN:	BC053039 , NM_010837 , NM_010837.1 , NM_010837.2 , NM_010837.3 , BE951502 , BE980868 , BQ444318 , BY259340
UniProt ID:	Q7TSJ2
Synonyms:	2810411E12Rik; Map-6; Mtap6; STOP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)