

Product datasheet for **MP208176**

Map7 Mouse qPCR Primer Pair (NM_008635)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	17761
Forward Sequence:	AGGTAGAGGAGGTCACAGTCGA
Reverse Sequence:	GGAACTACACCGACAGTCACAG
ACCN:	BC052637 , NM_008635 , NM_008635.1 , NM_008635.2 , BC042771 , BC044733 , BC046555 , BY277244
UniProt ID:	O88735
Synonyms:	E-MAP-115; MAP-7; mshi; mste; Mtap7; R75000; ste
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)