

Product datasheet for **MP208178**

Map9 Mouse qPCR Primer Pair (NM_001081230)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	213582
Forward Sequence:	ACAGTGAAGCCTACACCCAGGA
Reverse Sequence:	CTCCTCAGTGTGGCTGCTCTTT
ACCN:	<u>NM_001081230</u> , <u>NM_001081230.1</u> , <u>BC051460</u> , <u>BC061216</u> , <u>BC087935</u> , <u>BC148584</u> , <u>BC156817</u>
UniProt ID:	<u>Q3TRR0</u>
Synonyms:	5033421J10Rik; 5330427D05Rik; ASAP; Mtap9
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)