

Product datasheet for **MP208182**

Mapre3 Mouse qPCR Primer Pair (NM_133350)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	100732
Forward Sequence:	GCATGGGTCAATGACTCTCTGC
Reverse Sequence:	GCCTGGAACCTTGACCTTCCTCA
ACCN:	BC057918 , NM_133350 , NM_133350.1 , NM_133350.2 , BC048833 , BE948386 , BQ570892 , BY237685 , BY270434
UniProt ID:	Q6PER3
Synonyms:	AI790651; EB2; EB3; EBF3; RP3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)