

Product datasheet for **MP208195**

Mcm2 Mouse qPCR Primer Pair (NM_008564)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	17216
Forward Sequence:	CCGTTCCAAGGATGCCATTCTC
Reverse Sequence:	TGGAAGCCGTTGGCGGTGTTA
ACCN:	BC055318 , NM_008564 , NM_008564.1 , NM_008564.2 , BY324067
UniProt ID:	P97310
Synonyms:	AA959861; AW476101; BM28; CDCL1; Mcmd2; mKIAA0030
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)