

Product datasheet for **MP208212**

Slc25a53 Mouse qPCR Primer Pair (NM_001082412)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67062
Forward Sequence:	GCTTGCTGTACCTGCTGGGTTA
Reverse Sequence:	TTCGTTCTGGAGCTTCTGCTC
ACCN:	<u>NM_001082412</u> , <u>NM_001082412.1</u> , <u>NM_001082412.2</u> , <u>BC024604</u> , <u>BC152775</u>
Synonyms:	2310046F18Rik; 2810402A17Rik; AI448995; AI841128; Mcart6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)