

Product datasheet for **MP208234**

Mrpl24 Mouse qPCR Primer Pair (NM_026591)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67707
Forward Sequence:	CCATCTCTGATGAAGACTGGCAC
Reverse Sequence:	AGTTTCGCTGCCGAACGACTTG
ACCN:	BC004736 , BC025506 , BC031730 , NM_026591 , NM_026591.1 , NM_026591.2 , NM_026591.3
UniProt ID:	Q9CQ06
Synonyms:	2010005E08Rik; 2810470K06Rik; 6720473G22Rik; AA407670
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

