

Product datasheet for **MP208236**

Mrpl28 Mouse qPCR Primer Pair (NM_024227)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	68611
Forward Sequence:	GCTACAGATACGCCAACAACGAC
Reverse Sequence:	CACGGTGAACCTTCTGTCCAGG
ACCN:	<u>NM_024227</u> , <u>NM_024227.1</u> , <u>NM_024227.2</u> , <u>NM_024227.3</u> , <u>BC145680</u> , <u>BC002293</u> , <u>BC005797</u> , <u>BC080753</u> , <u>BC100322</u> , <u>BC145522</u> , <u>BC145682</u> , <u>BY138937</u>
UniProt ID:	<u>Q9D1B9</u>
Synonyms:	1110015G04Rik; L28mt; MAA; MAAT1; MRP-L28; p1; p15
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)