

Product datasheet for **MP208247**

Mrpl4 Mouse qPCR Primer Pair (NM_023167)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66163
Forward Sequence:	AGTCTCTCCGAGGCTTTGAGCA
Reverse Sequence:	CCTGAAGTTCCTCTGCCAGATG
ACCN:	<u>BC039983</u> , <u>BC061095</u> , <u>NM_023167</u> , <u>NM_023167.1</u> , <u>NM_023167.2</u> , <u>BC023078</u> , <u>BC037064</u> , <u>BC051685</u>
UniProt ID:	<u>Q9DCU6</u>
Synonyms:	1110017G11Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)