

Product datasheet for **MP208254**

Mrpl47 Mouse qPCR Primer Pair (NM_029017)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	74600
Forward Sequence:	GGTGCTTGGAGAAGGGACATCT
Reverse Sequence:	GCGATCCACATAAGGCATTGCG
ACCN:	BC087953 , NM_029017 , NM_029017.1 , NM_029017.2 , BC029173 , BC038065 , BY013327 , BY014313
UniProt ID:	Q8K2Y7
Synonyms:	4833424P18Rik; CGI-20; CGI-204; Gm9859; L47mt; MRP-L47; MTF/L47; NCM; NCM1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)