

Product datasheet for **MP208284**

Mrps31 Mouse qPCR Primer Pair (NM_020560)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57312
Forward Sequence:	AGGCGTTTGCTGATGAACCACC
Reverse Sequence:	TCTCAAAGGGCTGGTCAGCTAC
ACCN:	<u>BC018475</u> , <u>BC027430</u> , <u>NM_020560</u> , <u>NM_020560.1</u> , <u>NM_020560.2</u> , <u>BC006764</u> , <u>BY237441</u>
UniProt ID:	<u>Q61733</u>
Synonyms:	1500002D03Rik; AA960481; Imogen38; Imogen44
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)