

Product datasheet for **MP208285**

Mrps33 Mouse qPCR Primer Pair (NM_010270)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	14548
Forward Sequence:	GACTGGTATCCAAATCACAACACG
Reverse Sequence:	CCACGGAGCTTCTTTAGACGTC
ACCN:	BC092135 , BC096397 , NM_010270 , NM_010270.1 , NM_010270.2 , BC028645 , BC087881 , BU695761 , BY296250
UniProt ID:	Q9D2R8
Synonyms:	AI841153; CGI-139; Gdap3; MRP-S33; PTD003
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)