

Product datasheet for **MP208319**

Mapk8 Mouse qPCR Primer Pair (NM_016700)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26419
Forward Sequence:	CGCCTTATGTGGTACTCGCTA
Reverse Sequence:	TCCTGAAAGAGGATTTTGTGGC
ACCN:	<u>BC053027</u> , <u>NM_016700</u> , <u>NM_016700.1</u> , <u>NM_016700.2</u> , <u>NM_016700.3</u> , <u>NM_016700.4</u> , <u>BB619635</u>
UniProt ID:	<u>Q91Y86</u>
Synonyms:	AI849689; JNK; JNK1; Prkm8; SAPK1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)