

Product datasheet for **MP208323**

Mapk9 Mouse qPCR Primer Pair (NM_207692)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26420
Forward Sequence:	GTCAGTGGGTTGCATCATGGGA
Reverse Sequence:	ACTCTGCGGATGGTGTCCTAG
ACCN:	<u>BC028341</u> , <u>NM_207692</u> , <u>NM_207692.1</u> , <u>NM_207692.2</u> , <u>BB362351</u> , <u>BC024514</u>
UniProt ID:	<u>Q9WTU6</u>
Synonyms:	AI851083; JNK2; p54aSAPK; Prkm9
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)