

Product datasheet for **MP208338**

Map3k13 Mouse qPCR Primer Pair (NM_172821)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	71751
Forward Sequence:	GGAATGAGCCTGTCTCGGAGAA
Reverse Sequence:	CAAATGATGGCGGAAGAGTCCAC
ACCN:	<u>NM_172821</u> , <u>NM_172821.2</u> , <u>NM_172821.3</u> , <u>BB455895</u> , <u>BB549158</u> , <u>BP761466</u>
UniProt ID:	<u>Q1HKZ5</u>
Synonyms:	C130026N12Rik; LZK
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)