

## Product datasheet for **MP208346**

### Tab1 Mouse qPCR Primer Pair (NM\_025609)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66513
Forward Sequence:	TTCCTGGTGCTGATGTCAGAGG
Reverse Sequence:	AGGTCTGCTTGGCAAACCTCGGT
ACCN:	<a href="#">BC054369</a> , <a href="#">NM_025609</a> , <a href="#">NM_025609.1</a> , <a href="#">NM_025609.2</a> , <a href="#">BC027054</a> , <a href="#">BC041110</a> , <a href="#">BC114637</a> , <a href="#">BC134402</a> , <a href="#">NM_025609.3</a>
UniProt ID:	<a href="#">Q8CF89</a>
Synonyms:	2310012M03Rik; b2b449Clo; Map3k7ip1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)