

Product datasheet for **MP208354**

Map4k4 Mouse qPCR Primer Pair (NM_008696)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26921
Forward Sequence:	GCAGAATAGCCAAGCAGGACAG
Reverse Sequence:	TGAGTTGCTGGACCCTGAAGAG
ACCN:	<u>NM_008696</u> , <u>NM_008696.1</u> , <u>NM_008696.2</u> , <u>BB651050</u> , <u>BC056650</u> , <u>BC131673</u> , <u>BC141109</u> , <u>BC145390</u> , <u>BC145391</u> , <u>BC168387</u> , <u>BF721981</u> , <u>BQ571658</u>
UniProt ID:	<u>P97820</u>
Synonyms:	9430080K19Rik; AU043147; AU045934; AW046177; HGK; Nik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)